OE 7200

BAC A ULTRASONS 3L 100W 40KHz
ULTRASONIC CLEANER 3L 100W 40KHz

OE 7201

BAC ULTRASONS 9L 200W 40khz
ULTRASONIC CLEANER 9L 200W 40KHz

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
FEATURES

Normal and soft power control to resolve the cleaning blind spot .
 Degas function for increase the cleaning effect.
 Automatically ultrasonic frequency tracking.
 User-friendly and clear panel.
 Memory function.
 LED display for temperature and timer.
 High-performance transducers.
 Patent products.

Thank you for purchasing the CLAS Equipements ultrasonic cleaner. Please take the time to read these operating instructions before use and retain them for future reference. Failure to follow these instructions may lead to serious personal injury and damage to property.

INTRODUCTION

Ultrasonic cleaning is based on the cavitation effect caused by high frequency ultrasonic wave vibration signal in the fluid. Microscopic bubbles are formed, and they implode violently causing the cavitation which create an intense scrubbing action on the surface of the item being cleaned. The bubbles are small enough to penetrate microscopic crevices, cleaning them thoroughly and consistently. Ultrasonic cleaning is extremely effective at removing dirt and grime which would normally require tedious manual cleaning by hand. It has been used to clean a wide variety of instruments and mechanical parts such as carburetors, returning them to almost like new condition without damage to delicate parts.

PREPARATION

1. Carefully unpack the cleaner and remove all traces of packing materials from it. Visually inspect the cleaner for any parts that may have become loose or damaged during transit.

Contents :

A : Main machine B : Soundproof lid C : Power cord
 D : Outlet filter E : Mesh basket F : User manual

2. Place the cleaner on a flat, clean surface where the cooling fans will receive adequate ventilation and ensure all controls are set to off. Also ensure the drain tap is in the closed position.
3. Ensure the power cord is securely plugged into the cleaner and that no part is likely to come into contact with moisture.
4. Carefully fill at least $\frac{1}{2}$ of the tank with a solvent solution. For demanding cleaning we recommend a small amount of washing up liquid, this will help increase the cleaning performance. The cleaner is now ready for use.

ATTENTION

While the machine is working normally, ultrasonic and tank syntonny gives a well-proportioned sound, and no shudder on the surface of the water, yet there is spray made by the tiny bubbles. If there are discontinuous surges, please add or release a little of washing solution in the tank, stopping the surges is better for the objects cleaning.

SAFETY PRECAUTIONS

Keep awy from children !

Please read the following very carefully as failure to to comply may invalidate your guarantee.

TO AVOID ELECTRICAL SHOCK :

1. DO NOT run the cleaner continuously for more than one hour at a time, as doing so can damage the internal components.
2. DO NOT operate the unit without fluid in the tank. Always ensure the fluid is no higher than the max mark, Always ensure there is at least 1/2 solution.
3. DO NOT drop any item into the tank as this may cause damage to the transducer. Always place the items gently into the tank and use the basket whenever possible.
4. The more items that you place in your cleaning bath the less efficient it will clean. It is not advised to overlap items. Always allow plenty of clear space between the items.
5. Do keep the lid on during use. This will prevent splashes and reduce evaporation of the fluid.
6. Never immerse the machine or power cord in water or other liquid.
7. DO NOT touch the power plug with wet hands, especially when inserting or removing the plug.
8. DO NOT touch the unit if the machine has fallen into water during operation. Remove the power plug from the socket first.
9. DO NOT disassemble the machine, except by professionals.
10. UNPLUG the power source while filling or emptying the tank.
11. DO NOT spray water or liquid over the device and the control panel.
12. DO NOT operate the cleaner without proper grounding.
13. DO NOT place the device on a soft surface, where the vents could be blocked.
14. Upon completion of the cleaning cycle, turn the heater knob off and isolate the machine from the electrical supply.
15. Take care when adding or removing items from the cleaning tank as the fluid is likely to be hot and displaced fluid can damage the internal components. Any displaced fluid must be dried up immediately.

OE 7200

ULTRASONIC CLEANER 3L 100W 40KHz

OE 7201

ULTRASONIC CLEANER 9L 200W 40KHz

APPLICATION

Provided product is non-porous and can normally be immersed in water almost anything can be thoroughly cleaned.

Items Not Suitable For Ultrasonic Cleaning

Soft Jewelry: Pearls, emerald, ivory, coral. Agate, sea turtle shells, etc.

Welded, Plated and Glued Items: Welded or plated metal items, glued items.

Watches: Except divers watches with depth rating over 50m(150ft).

Others: Ceramic, camera filters with preexisting cracks.


UNIT FRONT VIEW/SIDE VIEW

- A : Tank
- B : Plastic carrying handles
- C : Draining
- D : Warning marking
- E : Soundproof lid
- F : Anti-slide feet


UNIT BACK VIEW

- A : Louver
- B : Socket
- C : Power switch
- D : QC PASSED sticker
- E : Rating label

CLAS®

FRONT PANEL

A : LED display temperature.

B : LED display cleaning time.

C : Heating button : Press this button for start/stop heating.

D : IncreaseTemperature : Press this button to set the temperature. Press one time to increase 1°C,holding the key to increase by 10°C

E : Temperature reduce: Press this button to set the temperature. Press one time to reduce 1°C, holding the key to reduce by 10°C.

F : Degas: Fresh solutions with many air which reduce effective ultrasonic action. Using Degas mode speeds up the degassing process and improves cleaning efficiency.

G : Soft : 60% power for delicate applications.

H : Normal 100% power to the tank for normal applications.

I : Timer increase : Press this button to set the cleaning time Press one time to increase 1minute,holding the key to increase by 10minutes.

J : Timer reduce : Press one time to reduce 1minute,holding the key to reduce by 10 minutes.

INSTRUCTIONS

Step	Action
1.	Select your cleaning solution, fill at least 1/2 of solution into the tank. Do not exceed the filling line. Note: Do not ever use alcohol, gasoline or flammable solutions. Doing so could cause a fire or explosion. Use only water-based solutions.
2.	Place the items into a basket, Slowly lower the basket into tank, do not allow items to contact the tank bottom (solution should cover the items).
3.	A : Plug the cleaner into grounded outlet. B : Power on and turn on the power switch on back side of the machine. Note: After the power on, the temperature you see on the machine is environment temperature. Time and temperature default as the last setting ones. (3 minutes default time for initial use)
TIMER 	Press and to setting the cleaning time. When timer count down to 00:00, the ultrasonic stop working. Note: The device with memory function and acquiesce in the time you setting last time.
HEATING 	Press to start/stop heating. Press and to setting the temperature. Note: 1. To get the best cleaning result, the solutions temperature we suggest is 40-60°C. Warm water and dish washing liquid can soften grease and improve cleaning result. 2. When setting the temperature, if the setting temperature exceed the environment temperature, press key, the heating working, and indicator on. If below environment temperature, the heating couldn't start, and indicator turn off. 3. Heating won't shut off automatically, if you don't press 4. The device with memory function and acquiesce in the temperature you setting last time..

INSTRUCTIONS

Étapes	Instructions
DEGAS MODE 	After timer and temperature setting, press , the corresponding indicator on and the device will work in degas mode : working 6 seconds, and stop 2 seconds, then working 6 seconds, cycle like this with countdown display. Intermittent operation of ultrasonic power for a quicker degassing of the cleaning liquid. If you want to transfer to other mode, just press  or , then the degas indicator turn off and indicator for normal or soft turn on. If want to stop working, press  again, then both heating and ultrasonic will stop and the corresponding indicator turn off .
NORMAL MODE 	After timer and temperature setting, Press  Normal key, the corresponding indicator on and the device will work in normal mode: continuously working with full ultrasonic power and countdown display. If you want to transfer to other mode, just press  or , the normal indicator turn off and indicator for degas or soft turn on; if you want to stop working, press  again, then both heating and ultrasonic will stop and the corresponding indicator turn off.
SOFT MODE 	After timer and temperature setting, press , the corresponding indicator on and the device will work in soft mode. If you want to transfer to other mode, just press  or , the soft indicator turn off and indicator of degas or normal turn on; if you want to stop working, press  again, then both heating and ultrasonic will stop and the corresponding indicator turn off.
LAST	After the cleaning process is finished,unplug the power cord,open the lid and remove the object. Finally, pour out the water after cooling. Clean and wipe the tank.

WARNING : Strong acid or alkaline cleaning solution will cause corrosion, rust and even puncture of tank or machine body. To overcome this problem, dilute to mild PH solution or request for tank made of specific grade of stainless steel .

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If you need components or parts, please contact the reseller.
In case of problems, please contact your authorized technician.